

New Generation Endoscope System





# FUJIFILM



FUJIFILM Corporation 26-30, NISHIAZABU 2-CHOME, MINATO-KU, TOKYO 106-8620, JAPAN http://www.fujifilm.com/products/medical/

SGE-157-00 Printed in JAPAN 201607-FP





# Welcome to **BLI&LCI World**





Linked Color Imaging

LMH	PUMP 👩 c	DN LIGHT	 POWER
● 1 2 3 HT LIMIT. LIGHT MO		LUMI.	0
<u>}</u> ↓		¢ ^	

This is a simulated image. Actual image will differ

# Superior Diagnosis Superior Usability









Expanding processing increases the color contrast nearby mucosal "Redness".

## 1 BLI (Blue Light Imaging)



# Colon

White Light mode





BLI mode





BLI mode

High contrast images suitable for observing microvascular and microsurface pattern are provided. Magnifying endoscopy is exellent with BLI.

White Light

Stomach



White Light mode



White Light mode

Bright, sharp, and stereoscopic images are accomplished with similar color tones to Xenon light source. Mega-pixel CMOS enables high-definition and quite low-level noise compatible.

Stomach



White Light mode

Esophagus



White Light mode

LCI would be helpful for detection with surface pattern and vessels. Slight color difference is visualized with natural tone, using "Red" component.

#### Observation modes can be switched by scope button.







LCI mode



LCI mode

Scope button "2" enables observation modes to be switched in the default setting.



BLI **BLI-bright** 



#### **3** Wide 170° field of view



## 4 High-definition images- Megapixel CMOS



Full HD display

## 5 Multi Zoom

#### Zoom function



\*When using a 26 inch HD LCD monitor

#### **Magnification Images**





Wide 170° field of view is available with EC-760R-V. Even areas that are hard to observe such as the reverse side of folds could be observed and approached smoothly.

High-definition images with quite low noise level are established by Megapixel CMOS sensor. It allows superior visualization for Full HD display.

```
n
```

Multi Zoom function are equipped on EG-760Z / EC-760ZP-V. With Continuous mode, Step Zoom mode of "2 Step", "3 Step" and "5 Step" are available. In this modes, images can be magnified in stages by simple press of button.

## **Superior Usability**



## **1** Variable Stiffness



Variable stiffness is equipped on EC-760R-V and EC-760ZP-V. The stiffness of the flexible portion can be adjusted with adjustment ring.

**Advanced Force Transmission** 

The flexible portion is designed to transmit operator's movements, pushing, pulling and rotating, to the distal end of endoscope.



Passing the sigmoid colon

#### **Adaptive Bending**

The end of flexible portion is soft, allowing the scope to bend with the angulations. Flexible portion is elastic, and easy to return to its straight shape after passing through angulations.



Passing the angulations





In deep insertion



After passing through angulations

#### **2** G7 control portion

G7 control portion is developed from ergonomics point of view.

Scope has a rounded surface to fit the hand, and button layout makes intuitive operation possible.





Identification color of instrument channel size Color of G7 control portion scopes nstrument channel diameter Model No. Corporate brand logo



Labels attached to the control portion, which show the model number and instrument channel inlet size.

## **4** Wide compatibility to conventional endoscope

Compatible with 700 series, 600 / 500 series endoscopes.



&

700 Series

	BLI	BLI-bright	LCI	WLI	FICE
700 Series	0	0	0	0	0
600 / 500 Series	×	×	×	0	0

#### **3** One Step Connector with Contact-free Technology



Scopes can be connected to light source in just 1 step operation. Scope cable connection is no longer required in setting up. One Step Connector enhances efficiency of clinical workflow.

#### **Contact-free Technology**

This's the generic name of below 3 points. It means connectors do not need to touch to transmit power and image data. By this technology, durability and reliability of scopes is expected to improve.





When compared to standard xenon light sources, the LED light source consumes about a third of the energy and lasts longer. Life time of the 4 LED light is expected for 6 years based on Fujifilm evaluation condition. Note: The warranty period is 1 year after date of purchase.





Conventional 600 / 500 Series



600/500 endoscopes can be used with White light and FICE mode. \* FICE : Flexible spectral Imaging Color Enhancement

#### Specification

#### Upper G.I. tract scopes

#### EG-760R



Field of view	140°
Observation range	2~100 mm
Bending capability	UP:210° DOWN:90° RIGHT:100° LEFT:100°
Working length	1,100 mm
Total length	1,400 mm
Distal end diameter	9.2 mm
Flexible portion diameter	9.3 mm
Minimum instrument channel diameter	2.8 mm

Product name: Video Endoscope GMDN: 38805 Generic name: Flexible video gastroduodenoscope



– Light guide

Light guide Objective lens

100° LEFT

nozzle

100° LEFT 100° RIGHT

#### EG-760Z



Field of view	Normal: 140° Close: 56°
Observation range	1.5~100 mm Normal: 3~100 mm Close: 1.5~2.5 mm
Bending capability	UP:210° DOWN:90° RIGHT:100° LEFT:100°
Working length	1,100 mm
Total length	1,400 mm
Distal end diameter	9.9 mm
Flexible portion diameter	9.8 mm
Minimum instrument channel diameter	2.8 mm

Normal: 140° Close: 56°

RIGHT: 160° LEFT: 160°

11.7 mm

11.8 mm

3.2 mm

1,330 mm (M) /1,690 mm (L)

1,650 mm (M) /2,010 mm (L)

1.5~100 mm

Normal: 3~100 mm Close: 1.5~2.5 mm UP: 180° DOWN: 180°

Product name: Video Endoscope GMDN: 38805 Generic name: Flexible video gastroduodenoscope

#### Lower G.I. tract scopes

#### EC-760R-V/M, I, L



#### EC-760ZP-V/M, L



Field of view	170°		
Observation range	2~100 mm		
Bending capability	UP: 180° DOWN: 180° RIGHT: 160° LEFT: 160°		
Working length	1,330 mm (M) /1,520 mm (l) / 1,690 mm (L)		
Total length	1,650 mm (M) /1,840 mm (l) / 2,010 mm (L)		
Distal end diameter	12.0 mm		
Flexible portion diameter	12.0 mm		
Minimum instrument channel diameter	3.8 mm		
Product name: Video Endoscope			

GMDN: 36117 Generic name: Flexible video colonoscope

Field of view

Observation range

Bending capability

Distal end diameter

Minimum instrument channel diameter

Flexible portion diameter

Working length

Total length





Light guide Objective len



VP-7000

Digital Output	DVI (Resolution: 1280 × 1024 pixels, 1920 × 1080 pixels) HD-SDI (Resolution: 1920×1080 pixels)			
Output Connector	DVI-D: DVI-I: HD-SDI:	2 channel 1 channel 2 channel	S VIDEO: VIDEO: RGB TV:	1 channe 1 channe 1 channe
Input Connector	1 channel PoP			
Control Connector	Light Source I/F (37p): Light Source I/F (Mini D-sub 15p) REMOTE (BNC): Peripherals (D-sub 9p): Keyboard:	1 channel : 1 channel 2 channel 2 channel 1 channel	Card reader: Digital printer Footswitch: Network:	1 channe : 1 channe 1 channe 1 channe
Type of color	NTSC/PAL			
Iris	Auto / Peak / Ave			
Applicable Endoscope	700/600/500 series endosc	opes		
Power Rating	Voltage: 100 Frequency: 50/6 Current consumption: 0.8-0	to 240 V 60 Hz 0.5 A		
Dimensions( $W \times H \times D$ )	390 × 110 × 485 mm (inclu	ding projectio	ר)	
Weight	9.0 kg			

Product name: Processor GMDN: 18034 Generic name: Endoscopic Video image processor

VP-7000

BL-7000

#### New Accessories (Valve, Tank)

#### For routine examination



-

Air / Water Valve AW-603

Suction Valve SB-605

#### Used with CO<sub>2</sub> Regulator "GW-100"



Product name:Video Endoscope GMDN: 36117 Generic name: Flexible video colonoscope

#### BL-7000

Light Source	LED Maximum light output: 1400 lm		
Light Control	Autmatic light control by the control signal from the video processor		
Light Cooling Method	Forced air cooling		
Air Supply Pump	Available at 4-levels (Hi / Mid / Low / Off)		
Power Rating	Voltage: 100 to 240 V ~ Frequency: 50 / 60 Hz Current consumption: 1.2-0.7A		
Dimensions ( $W \times H \times D$ )	$390 \times 155 \times 485$ mm (including projection)		
Weight	12 kg		

Product name: Light source GMDN: 35158 Generic name: Endoscopic light source, line-powered





Water Tank WT-603



Water Tank WT-604G