



BioVac® direct suction device...

Specializing in evacuation and lavage of the GI tract, the BioVac® device is designed to suction hard-to-remove, viscous materials, like gelatinous blood and stringy clots, which may otherwise clog the control head of the endoscope. By bypassing the scope control head, the BioVac® device alleviates the intraprocedural clogging of the scope suction valve, while offering considerable suction

The BioVac® device offers...

power and volume.

- Enhanced suction, facilitating increased volume* and quicker, more powerful evacuation.
- The ability to irrigate as much as desired instilling large volumes of fluid is no longer a nuisance with the BioVac® device's quick evacuation.
- Direct visualization to easily target cleansing and evacuation sites, allowing clinicians to focus on the patient's diagnosis and/or treatment.
- Instrument access, enabling the clinician to utilize the BioVac® device's suction capabilities as well as provide therapeutic treatment via the accessory channel (see table below for devices featuring instrument access).

Experience the power of the BioVac® device on your next GI bleed, colonic decompression, poor prep or retained residual food case.

trumpet valve for suction control

For more information on the BioVac® device, contact US Endoscopy today.

product number	description	scope compatibility	device access	y-port	units/box
00711511	direct suction device	Pentax	yes	no	5
00711512	direct suction device	Olympus/Fujinon**	yes	no	5
00711513	direct suction device	Olympus/Fujinon**	yes	yes	5

*The size of particles and volume of liquid suctioned is limited to the endoscope's channel size.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™

are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy Fujinon® is a registered trademark of Fuji Film Co., Ltd.

Olympus® is a registered trademark of Olympus Optical, Inc

PENTAX® is a registered trademark of Asahi Kogaku Kogyo Kabushiki Kaisha



listening...and delivering solutions®

5976 Heisley Road Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopv.com

760410 Rev C