

A large, detailed photograph of an underwater scene. In the foreground, a large sea turtle swims towards the left. In the background, another sea turtle swims away. The water is clear blue, and sunlight filters down from the surface, creating a bright, shimmering effect. A coral reef is visible at the bottom of the frame.

medwork®
for mindful people

Biopsy







Product overview

www.medwork.com

Biopsy

EndoBite² Biopsy forceps



	Product Code	Order No.	Characteristics	ø Jaws	Length	ø Working Channel	Unit
	BI01-C4-18-120	502004	Oval jaws, coated, Bronchoscopy	1,8 mm	120 cm	2,0 mm	10
	BI01-C4-18-180	502005	Oval jaws, coated, Paediatrics	1,8 mm	180 cm	2,0 mm	10
	BI01-C4-18-260	502006	Oval jaws, coated, Balloon Enteroscopy	1,8 mm	260 cm	2,0 mm	10
	BI01-C4-23-180	501900	Oval jaws, mono pull wire, coated	2,3 mm	180 cm	2,8 mm	10
	BI01-C4-23-230	501901	Oval jaws, mono pull wire, coated	2,3 mm	230 cm	2,8 mm	10
	BI01-C4-30-230*	502558	Oval jaws, mono pull wire, coated	3,0 mm	230 cm	3,2 mm	10
	BI01-D4-23-180	501898	Oval jaws, mono pull wire, coated, spike	2,3 mm	180 cm	2,8 mm	10
	BI01-D4-23-230	501899	Oval jaws, mono pull wire, coated, spike	2,3 mm	230 cm	2,8 mm	10
	BI01-D4-30-230*	502241	Oval jaws, mono pull wire, coated, spike	3,0 mm	230 cm	3,2 mm	10
	BI01-C3-23-180	500160	Oval jaws, mono pull wire	2,3 mm	180 cm	2,8 mm	10
	BI01-C3-23-230	500161	Oval jaws, mono pull wire	2,3 mm	230 cm	2,8 mm	10
	BI01-D3-23-180	500164	Oval jaws, mono pull wire, spike	2,3 mm	180 cm	2,8 mm	10
	BI01-D3-23-230	500165	Oval jaws, mono pull wire, spike	2,3 mm	230 cm	2,8 mm	10
	BI01-G4-23-180	502162	Alligator, mono pull wire	2,3 mm	180 cm	2,8 mm	10
	BI01-G4-23-230	502163	Alligator, mono pull wire	2,3 mm	230 cm	2,8 mm	10
	BI01-H4-23-180	502556	Alligator, mono pull wire, spike	2,3 mm	180 cm	2,8 mm	10
	BI01-H4-23-230	502557	Alligator, mono pull wire, spike	2,3 mm	230 cm	2,8 mm	10

*Available from May 2019

EndoBite² with mono pull wire

Ensures large sizes of sample material ideal for histological assessment because the cut surfaces are completely closed, even if the position of the endoscope is challenging. The mono pull wire prevents the branches from tilting and thus protects the working channel against wear and damage.



Biopsy

NEW
WITH COLOUR CODE

Gastro Colo Broncho DBE



Spike No spike

ø 1,8 mm ø 2,3 mm ø 3,0 mm



With medwork BIOpin
Easy removal of biopsies from the jaws



contact@medwork.com
www.medwork.com

Tel.: +49 (0) 91 93 / 500 900
Fax: +49 (0) 91 93 / 500 999

medwork GmbH
Medworkring 1
91315 Höchstadt/Aisch
Germany