



Instruments
and accessories

for vessel sealing

Legend

Important notes


The authors and Erbe Elektromedizin GmbH have taken the greatest possible care while compiling the catalogue with electrosurgical accessories. Nevertheless, it is impossible to exclude the possibility that errors may be contained herein.

The information and recommendations given here may not be construed as constituting any basis for any claims against Erbe Elektromedizin GmbH. Should legal regulations stipulate liability, such liability shall be limited to intentional misconduct or gross negligence.

The data on the use of instruments is based on clinical experience. Medical science is subject to a process of permanent change based on research and clinical experience. Therefore the products are subject to alterations in design and technique. Because of production tolerances dimensions in the catalogue may marginally differ.

We confirm that none of our products contain latex, except the rubber straps, Art. No. 20592-009 and 20592-011.

Further product information



Please read carefully the corresponding notes on use before using any instrument.

Vessel sealing

Instruments and accessories

- 04 BiClamp
- 10 BiCision
- 11 BiClamp LAP forceps
- 13 BiClamp E LAP forceps
- 15 Accessories for treatment for BiClamp LAP and BiClamp E LAP forceps

Overview

Field of application

Open surgery

Laparoscopy

Open surgery and laparoscopy

Shape of instrument


BiClamp

BiClamp LAP and BiClamp E LAP

BiCision

Electrosurgical unit

VIO 3



BiClamp



BiClamp



210 mm

260 mm

BiClamp



270 mm

280 mm

100 %







Instruments and accessories

BiClamp

BiClamp



No.	
= 1	0.25 kVp
	max 95 °C
30195-122	max 138 °C

Product type	Length	Shape	Jaws	No.
BiClamp 150 C 	150 mm	bent 23°	smooth	20195-221
BiClamp 201 T 	200 mm	bent 18°	smooth	20195-202
BiClamp 210 	210 mm	bent 25°	corrugated	20195-200
BiClamp 260 C 	260 mm	bent 18°	smooth	20195-299
BiClamp 271 T 	270 mm	bent 18°	smooth	20195-203
BiClamp 280 	280 mm	bent 25°	smooth	20195-280

100 %

Instruments and accessories

BiClamp

BiClamp 150 C, bent 23°, smooth, length 150 mm

with connecting cable 4 m and MF plug, with ceramic coating, e.g. for thyroidectomy



No. 20195-221



= 1

0.25 kVp

max 95 °C



30195-122

max 138 °C

BiClamp 201 T, bent 18°, smooth, length 200 mm

with connecting cable 4 m and MF plug, with thermal insulation, for open surgical procedures, e.g. intestinal surgery



No. 20195-202



= 1

0.25 kVp

max 95 °C



30195-122

max 138 °C

BiClamp 210, bent 25°, smooth, length 210 mm

with connecting cable 4 m and MF plug, with corrugated surfaces, e.g. for vaginal hysterectomy



No. 20195-200



= 1

0.25 kVp

max 95 °C



30195-122

max 138 °C

Instruments and accessories

BiClamp

BiClamp 260 C, bent 18°, smooth, length 260 mm

with connecting cable 4 m and MF plug, with ceramic coating, e.g. for prostatectomy or pediatric surgery



No. 20195-299



= 1

0.25 kVp

max 95 °C



30195-122

max 138 °C

BiClamp 271 T, bent 18°, smooth, length 270 mm

with connecting cable 4 m and MF plug, with thermal insulation, for open surgical procedures, e.g. intestinal surgery



No. 20195-203



= 1

0.25 kVp

max 95 °C



30195-122

max 138 °C

BiClamp 280, bent 25°, smooth, length 280 mm

with connecting cable 4 m and MF plug. e.g. for open cystectomy or prostatectomy



No. 20195-280



= 1

0.25 kVp

max 95 °C



30195-122

max 138 °C

Instruments and accessories

BiCision

BiCision S, ø 5 mm, shaft length 200 mm

for open and laparoscopic procedures



No. 20195-310

= 5

0.25 kVp

30195-331



STERILE EO

BiCision M, ø 5 mm, shaft length 350 mm

for open and laparoscopic procedures



No. 20195-311

= 5

0.25 kVp


30195-331



STERILE EO

BiCision L, ø 5 mm, shaft length 450 mm

for open and laparoscopic procedures



No. 20195-312

= 5

0.25 kVp

30195-331



STERILE EO

Instruments and accessories

BiClamp LAP forceps

BiClamp LAP forceps, Maryland, semi-deep ribbed, shaft \varnothing 5 mm, non-adhesive coating, length 340 mm
with connecting cable 4 m and MF plug, complete instrument



No. 20195-134

= 1	0.25 kVp
	max 95 °C
30195-119	max 138 °C

BiClamp LAP forceps, Maryland, smooth, shaft \varnothing 5 mm, non-adhesive coating, length 340 mm
with connecting cable 4 m and MF plug, complete instrument



No. 20195-135

= 1	0.25 kVp
	max 95 °C
30195-119	max 138 °C

BiClamp LAP forceps, fenestrated, semi-deep, shaft \varnothing 5 mm, non-adhesive coating, length 340 mm
with connecting cable 4 m and MF plug, complete instrument



No. 20195-136

= 1	0.25 kVp
	max 95 °C
30195-119	max 138 °C

BiClamp LAP forceps, fenestrated, smooth, shaft \varnothing 5 mm, non-adhesive coating, length 340 mm
with connecting cable 4 m and MF plug, complete instrument



No. 20195-137

= 1	0.25 kVp
	max 95 °C
30195-119	max 138 °C

Instruments and accessories

BiClamp LAP forceps

BiClamp LAP forceps, Kelly, semi-deep, shaft \varnothing 5 mm, non-adhesive coating, length 340 mm

with connecting cable 4 m and MF plug



No. 20195-228



= 1

0.25 kVp

max 95 °C



30195-119

max 138 °C

BiClamp 110, smooth, shaft \varnothing 5 mm, non-adhesive coating, length 110 mm

with connecting cable 4 m and MF-plug, complete instrument,
e.g. for minimally invasive thyroidectomy (MIVAT)



No. 20195-230



= 1

0.25 kVp

max 95 °C



30195-096

max 138 °C

Single components		Piece	No.
Insert for BiClamp LAP forceps, Kelly, shaft \varnothing 5 mm, length 340 mm (for 20195-228)		1	20195-236
Insert for BiClamp LAP forceps, Maryland, semi-deep ribbed, shaft \varnothing 5 mm, length 340 mm (for 20195-134)		1	20195-146
Insert for BiClamp LAP forceps, Maryland, smooth, shaft \varnothing 5 mm, length 340 mm (for 20195-135)		1	20195-147
Insert for BiClamp LAP forceps, fenestrated, semi-deep, shaft \varnothing 5 mm, length 340 mm (for 20195-136)		1	20195-148
Insert for BiClamp LAP forceps, fenestrated, smooth, shaft \varnothing 5 mm, length 340 mm (for 20195-137)		1	20195-149
insert for BiClamp 110, smooth, shaft \varnothing 5 mm, length 110 mm (for 20195-230)		1	20195-232
Shaft tube for bipolar LAP forceps and BiClamp LAP forceps		1	20195-141
Shaft tube for BiClamp 110, \varnothing 5 mm, length 110 mm (for 20195-230)		1	20195-234
Handle for BiClamp LAP forceps and BiClamp 110, with connecting cable 4 m and MF plug		1	20195-145

Instruments and accessories

BiClamp E LAP forceps

BiClamp E LAP forceps, Maryland, semi-deep, shaft \varnothing 5 mm, length 340 mm

with connecting cable 4 m and MF plug



No. 20195-246



= 1

0.25 kVp

max 95 °C



30195-119

max 138 °C

BiClamp E LAP forceps, Maryland, smooth, shaft \varnothing 5 mm, length 340 mm

with connecting cable 4 m and MF plug



No. 20195-247



= 1

0.25 kVp

max 95 °C



30195-119

max 138 °C

BiClamp E LAP forceps, fenestrated, semi-deep, \varnothing 5 mm, length 340 mm

with connecting cable 4 m and MF plug



No. 20195-248



= 1

0.25 kVp

max 95 °C



30195-096

max 138 °C

BiClamp E LAP forceps, fenestrated, smooth, shaft \varnothing 5 mm, length 340 mm

with connecting cable 4 m and MF plug



No. 20195-249



= 1

0.25 kVp

max 95 °C



30195-096

max 138 °C

Instruments and accessories

BiClamp E LAP forceps

BiClamp E LAP forceps, Kelly, semi-deep, shaft \varnothing 5 mm, length 340 mm

with connecting cable 4 m and MF plug










No. 20195-229

0.25 kVp

max 95 °C

i 30195-119

max 138 °C

Single components		Piece	No.
Insert for BiClamp E LAP forceps, Kelly, shaft \varnothing 5 mm, length 340 mm (for 20195-229)		1	20195-236
Insert for BiClamp E LAP forceps, Maryland, semi-deep ribbed, shaft \varnothing 5 mm, length 340 mm (for 20195-246)		1	20195-146
Insert for BiClamp E LAP forceps, Maryland, smooth, shaft \varnothing 5 mm, length 340 mm (for 20195-247)		1	20195-147
Insert for BiClamp E LAP forceps, fenestrated, semi-deep, shaft \varnothing 5 mm, length 340 mm (for 20195-248)		1	20195-148
Insert for BiClamp E LAP forceps, fenestrated, smooth, shaft \varnothing 5 mm, length 340 mm (for 20195-249)		1	20195-149
Shaft tube for BiClamp E LAP forceps, \varnothing 5 mm, length 340 mm		1	20195-241
Handle for BiClamp E LAP forceps, with connecting cable 4 m and MF plug		1	20195-240

Instruments and accessories


Accessories for treatment for BiClamp LAP and BiClamp E LAP forceps


Cleaning brush for suction channel max. \varnothing 3.5 mm, length 650 mm



No. 20191-279

 = 1

 max 95 °C

 max 134 °C

Irrigation tube for shaft tube

Lap instruments



No. 20195-201

 = 1

Erbe Elektromedizin GmbH
Waldhoernlestrasse 17
72072 Tuebingen
Germany

Phone +49 7071 755-0
Fax +49 7071 755-179
info@erbe-med.com
erbe-med.com