



enteroscopy devices

Enteroscopy devices

US Endoscopy offers a full line of single-use, enteroscopy devices.



Reach® overtube

00712140

scope OD 10.5-11.7mm

length 80cm

overtube OD 16.5mm

overtube ID 13.5mm

The Reach® overtube helps splint the endoscope from the lower esophageal sphincter to the pylorus, offering a more direct path to reach proximal small bowel findings.

minimizes reformation of gastric loops

allows for further advancement of the scope

promotes one-to-one movement during advanced therapy



Overtube minimizes reformation of the loop

Carr-Locke injection needle

00711811

sheath diameter 2.5mm
length 230cm
needle projection 5mm
needle gauge 25

00711812

sheath diameter 2.5mm
length 230cm
needle projection 5mm
needle gauge 23

00711813 (Pentax)

sheath diameter 2.5mm
length 230cm
needle projection 5mm
needle gauge 25

00711814

sheath diameter 2.5mm
length 230cm
needle projection 4mm
needle gauge 23



The Carr-Locke needle is an ideal choice for injection therapy during enteroscopy procedures.

patented design ensures consistent needle projection every time in the most tortuous conditions

luer-lock, spring-loaded handle automatically retracts the needle, reducing the risk of inadvertent needle sticks and scope damage

Articulator™ injection needle

00711808

sheath diameter 2.5mm
length 350cm
needle projection 5mm
needle gauge 25

The Articulator™ needle features a kink resistant, stainless steel sheath to pass smoothly and perform in fully articulated endoscopes.



AdvanCE[®] capsule delivery device

00711144

sheath diameter 2.5mm

length 180cm

Allows for direct endoscopic placement of video capsules* in patients who are unable to swallow or pass the capsule through the pylorus.

facilitates endoscopic delivery of video capsules into the stomach or duodenum

enables small intestinal surveillance in challenging patients with dysphagia, pill phobia, gastroparesis, anatomical abnormalities, and more



*For use with the PillCam SB video capsule

Roth Net[®] retriever

00711152

sheath diameter 2.5mm

length 350cm

approximate net size 3x6cm

Designed for removal of excised polyps and foreign bodies within the convolutions of the small bowel.

engineered with a strong, reinforced net for extra security

transparent net material allows for continued examination with easy visibility



Cytology brush

0071**1405**

sheath diameter 2.4mm

length 350cm

bristle diameter 3.0mm

The unique bristle configuration of the enteroscopy length cytology brush helps the clinician obtain an abundant cell sample from the small bowel.



Biopsy forceps

standard

0071**1291**

sheath diameter 2.3mm

length 350cm

oval with spike

hot

0071**1295**

sheath diameter 2.5mm

length 350cm

oval

State-of-the-art manufacturing processes ensure that these biopsy forceps are sharp and well aligned, allowing for precise tissue samples every time.



Short Throw™ polypectomy snare

00711101

sheath diameter 2.4mm

length 350cm

approximate snare size 2.5x5.4cm

The Short Throw™ snare offers superior loop memory and the comfort of one-handed deployment.

short throw handle offers a 25% stroke travel reduction

designed to work in tortuous conditions

Bipolar hemostasis probes

00711841

size 7FR

length 350cm

dual plug

00711843

size 7FR

length 350cm

single plug

00711847

size 10FR

length 350cm

single plug

The hemostasis probe features a unique, Hydromer tip to help prevent tissue from sticking to the probe.

Jumbo channel cleaning brush

0071**1606**

sheath diameter 1.7mm
length 350cm
bristle diameter 7-9mm taper

Designed with a textured, non-slip catheter and soft, specially designed nylon bristles for maximum cleaning effort in scopes that range from 2.8mm and larger.



Jumbo channel cleaning brush

0071**1610** / 0071**1611**

bristle head length 4cm
bristle diameter 10mm

A full head of soft bristles to aid in removal of gross material, the ridged handle is designed with grooves for a comfortable, non-slip grip.



For more devices available for single balloon and double balloon enteroscopes (230cm length), please visit www.usendoscopy.com.



listening... and delivering solutions®

5976 Heisley Road
Mentor, OH 44060
phone (1) 440 / 639-4494
fax (1) 440 / 639-4495
customer service (1) 800 / 769-8226
email global@usendoscopy.com

www.usendoscopy.com